

REF 09000-110-E

Mobile Collection System

The MCS[®]+ 9000 Mobile Collection System is dedicated to blood component collection, including platelets and plasma.

It also performs therapeutic plasma exchange.



Technical Information

Machine (closed)	Depth: 13.38 inches Height: 17.25 inches Width: 22.25 inches
Machine (in operation)	Depth: 22.25 inches Height: 27 inches Width: 22.25 inches
Weight of machine	56 lbs
Construction	Sheet Metal
Electrical voltage	100/220 VAC +/- 10%
Operating frequency	47 - 63 Hz
Fuse rating	5.0 A domestic
Ambient operating temperature	+10° C to +27° C
Storage humidity level	8% to 80%, non-condensing

Accessories

Device Stand/Cart	S10219-00
Cuff Adult	52304-02
Cuff Large	52304-03
Mobile Transport Case	09012-00

Associated Products

SDPE Protocol portfolio	09081-00
TPE Protocol portfolio	09082-00

Ordering Information

09000-110-E	MCS+, 110V ENGL, LN 09000-110-E
-------------	---------------------------------



For a list of worldwide office locations and contact information, visit www.haemonetics.com/officelocations