


# Sensor-Guided TAVI




## SavvyWire®

Pre-Shaped Pressure Guidewire

**TAVI. Efficiency. Haemodynamics.**

The SavvyWire guidewire, the first Sensor-Guided TAVI solution, is designed to optimize TAVI through efficient, predictable wire performance, continuous haemodynamic measurement and LV pacing capabilities.



### TAVI

**High performance TAVI wire.** The SavvyWire guidewire is engineered for workhorse guidewire performance to support stable valve delivery and positioning.



### Efficiency

**Optimise without compromise.** The SavvyWire guidewire offers several advantages to streamline the procedure by reducing procedural steps and supporting case efficiency.



### Haemodynamics

**Continuous, invasive haemodynamic feedback.** Powered by Fidelis® technology, The SavvyWire guidewire delivers continuous, accurate haemodynamic measurement that enables efficiencies and support procedural decision-making for lifetime TAVI management.


**HAEMONETICS®**

Interventional Technologies

# SavvyWire®

## Pre-Shaped Pressure Guidewire

Designed to optimise the TAVI procedure



### Fidela® Fiber Optic Pressure Sensor Specifications

|                            |   |
|----------------------------|---|
| Operating pressure         | -30 to +300 mmHg  |
| Accuracy                   | ± 1 mmHg plus ± 1% of reading (-30 to 50 mmHg),<br>3% of reading (50 to 300 mmHg) |
| Zero thermal effect        | 0.3 mmHg / °C   |
| Sensitivity thermal effect | 0.3% / °C   |
| Zero drift                 | < 1 mmHg / hr   |

### Savvywire Pre-Shaped Pressure Guidewire Ordering Information

| Product Ref Number | Description       | Tip Length | Tip Diameter |
|--------------------|-------------------|------------|--------------|
| F3001              | SavvyWire X-Small | 2.9 cm     | 3.2 cm       |
| F3002              | SavvyWire Small   | 3.6 cm     | 4.2 cm       |

### OptoMonitor® 3 System Ordering Information

| Product Ref Number | Description                    |
|--------------------|--------------------------------|
| F2321              | OptoMonitor 3 Display Unit 15" |
| F2310              | OptoMonitor 3 Optical Unit     |
| F2330              | OptoMonitor 3 Handle Unit      |
| F232030            | Activation Key TAVI 1          |

Please consult product labels and instructions for use for indications, contraindications, warnings, precautions and adverse events. See SavvyWire IFU LBL-2015-03 and OptoMonitor IFU LBL-2019-41.

For a list of worldwide office locations and contact information, visit [www.haemonetics.com/contact-support](http://www.haemonetics.com/contact-support)

Learn More: [www.haemonetics.com](http://www.haemonetics.com)