

NEW

Simply Bigger



**With 58% More
Collagen***

VASCADE MVP[®] XL

VASCADE MVP® XL System

VASCADE MVP® and VASCADE MVP® XL Dry, Uncompressed Collagen Comparisons

With 58% More Collagen*




VASCADE MVP®



VASCADE MVP® XL



VASCADE MVP®



VASCADE MVP® XL

Features	VASCADE MVP® System	VASCADE MVP® XL System	Increase
Collagen	12mg ±3mg	19mg ±3mg	58%
Disc Size	23F	25F	9%
Sheath Size Compatibility	6-12F ID (15F OD)	10-12F ID (15F OD)	None

* As compared to VASCADE MVP

Please consult product labels, and instructions for use for indications, contraindications, warnings, precautions and adverse events.
See VASCADE MVP XL IFU 5611 at: info.haemonetics.com/xl
For a list of worldwide office locations and contact information, visit www.haemonetics.com/office/locations

© 2024 Haemonetics Corporation. Haemonetics, and VASCADE MVP are trademarks or registered trademarks of Haemonetics Corporation in the U.S., other countries or both. 05.2024 USA.
LIT 2072 Rev AA

HAEMONETICS®