

# Global Haemostasis

## Heparin Neutralisation Cartridge

Supporting key steps in protocol-based haemostasis management

**The Global Haemostasis–Heparin Neutralisation (GH-HN) cartridge aids clinicians in making timely protocol-driven decisions across the perioperative pathway.**

In addition to the standard Global Haemostasis assays, this cartridge offers the additional sensitivity to heparinisation for both the RapidTEG™ and Functional Fibrinogen assays.



### A comprehensive view of the GH-HN cartridge

#### What is the added value?

By reliably assessing coagulation in heparinised patients, GH-HN supports:

- The drive for a more effective programmatic based adoption of Viscoelastic Testing (VET)
- Blood product optimisation and utilisation
- Targeted and personalised transfusion options
- Simplified view of data for critical decisions making



#### Setting a new standard in Viscoelastic Testing

Together with the Global Haemostasis and the **PlateletMapping®** cartridges, the **GH-HN** cartridge supports the **TEG® 6s** analyser and **TEG Manager®** software to deliver a unique single platform that combines global haemostasis and platelet function with the most comprehensive clinical interpretation.

# Features of the TEG® 6s Viscoelastic Testing Platform

## Feature

Simple, self-contained cartridge system	■
Test preparation without complex pipetting and mixing	■
Smaller blood volume per test	■
Actionable results available within minutes	■
Comprehensive hemostasis analysis including PlateletMapping® assay	■
Measures same viscoelastic properties as previous generations of TEG® analyzers	■
Provides more information than routine coagulation tests	■
Smallest footprint and lightest weight viscoelastic testing platform on the market	■
Vibration resistant	■
Robust self-monitoring quality control	■
Device can be used in the lab or near the patient as indicated	■
Remote results viewing and device administration with TEG Manager® software	■
TEG Manager Interpretation Guidance helps simplify and optimize consistency of care	■



## TEG 6s platform tests available

Cartridge	Assays	Specimen Type
<b>Global Hemostasis</b>	Kaolin	Citrated whole blood
	Kaolin with Heparinase	
	RapidTEG™	
	Functional Fibrinogen	
<b>Global Hemostasis – HN</b>	Kaolin	Citrated whole blood
	Kaolin with Heparinase	
	RapidTEG™ with Heparinase	
	Functional Fibrinogen with Heparinase	
<b>PlateletMapping®</b>	PlateletMapping ADP & AA	Heparinized whole blood
	PlateletMapping ADP	



## Ordering Information

Description	Item Code	Description	Item Code
TEG 6s analyzer	07-690	PlateletMapping ADP & AA Cartridge (box/10)	07-614
Printer, USB	118034-00	PlateletMapping ADP Cartridge (box/10)	07-615
Printer Paper, USB (roll)	118035-00	Level 1 QC Test Kit (box/12)	07-664
Barcode Scanner	118332-00	Level 2 QC Test Kit (box/12)	07-665
Global Hemostasis Cartridge (box/10)	07-601		
Global Hemostasis-HN Cartridge (box/10)	07-604		

The products may not be available in all countries. To check availability in your region, please contact us for more information.

For a list of worldwide office locations and contact information, visit [www.haemonetics.com/officelocations](http://www.haemonetics.com/officelocations)